

ChiralTek Column Compatible and Alternative Replacement List

ChiralTek	Daicel	Supelco/Astec	Phenomenex	Merck/Astec
ChiralAM-2	IA		Lux i-Amylose-1	
ChiralCE-2	IB			
ChiralCE-4	IC		Lux i-Cellulose-5	
ChiralAM-13	ID			
ChiralAM-4	IE			
ChiralAM-3	IF			
ChiralAM-5	IG			
ChiralAM-6	IH			
ChiralCE-7	IJ			
ChiralCE-5	IK			
ChiralAM-1				
ChiralAM-7				
ChiralAM-8				
ChiralAM-9				
ChiralAM-10				
ChiralAM-11				
ChiralAM-12				
ChiralCE-1				
ChiralCE-3	IM			
ChiralCE-6				
ChiralCE-8				
ChiralCE-9				
ChiralCE-10				
ChiralCE-11				
ChiralCE-12				
ChiralCE-13				
ChiralCTCE-1	OC/OC-H			
ChiralCTCE-2	OD/OD-H		Lux Cellulose-1	
ChiralCTCE-7	OJ/OJ-H		Lux Cellulose-3	
ChiralCTCE-3	OZ/OZ-H		Lux Cellulose-2	
ChiralCTCE-9	OX/OX-H		Lux Cellulose-4	
ChiralCTCE-8	OG/OG-H			
ChiralCTAM-6	AS/AS-H			
ChiralCTAM-2	AD/AD-H		Lux Amylose-1	
ChiralCTAM-3	AZ/AZ-H			
ChiralCTAM-10	AY/AY-H		Lux Amylose-2	
ChiralCTCE-2R	OD-R/OD-RH			
ChiralCTCE-7R	OJ-R/OJ-RH			
ChiralCTCE-3R	OZ-R/OZ-RH			
ChiralCTCE-9R	OX-R/OX-RH			
ChiralCTAM-6R	AS-R/AS-RH			
ChiralCTAM-2R	AD-R/AD-RH			
ChiralCTAM-3R	AZ-R/AZ-RH			
ChiralCTAM-10R	AY-R/AY-RH			
ChiralCTCE-1R	OC-R/OC-RH			
ChiralCTCE-8R	OG-R/OG-RH			
ChiralCF6	CR(+)			
ChiralCF6	CR(-)			
Chiral β CD-0		Astec Cyclobond I 2000		
Chiral γ CD-0		Astec Cyclobond II		
Chiral β CD-2		Astec Cyclobond I 2000 DMP		
Chiral β CD-1/ β CD-2/ β CD-3/ β CD-4/ β CD-5/ β CD-6/ β CD-7/ β CD-8/ β CD-9/ β CD-10/ β CD-11/ β CD-12/ β CD-13				ChiralDex series LC columns

ChiralTek Column Compatible and Alternative Replacement List

ChiralTek	Daicel	Supelco/Astec	Phenomenex	Merck/Astec
Chiral α CD-0/ α CD-1/ α CD-2/ α CD-3/ α CD-4/ α CD-5/ α CD-6/ α CD-7/ α CD-8/ α CD-9/ α CD-10/ α CD-11/ α CD-12/ α CD-13				
Chiral γ CD-1/ γ CD-2/ γ CD-3/ γ CD-4/ γ CD-5/ γ CD-6/ γ CD-7/ γ CD-8/ γ CD-9/ γ CD-10/ γ CD-11/ γ CD-12/ γ CD-13				
The complex chiral selector-immobilized columns & media				
ChiralAMCE-1				
ChiralAMCE-2	IA + IB			
ChiralAMCE-3	IF +			
ChiralAMCE-4	IC + IE			
ChiralAMCE-5	IG + IK			
ChiralAMCE-6	IH +			
ChiralAMCE-7	IJ +			
ChiralAMCE-13	ID +			
...				
ChiralAM2AM3	IA + IF			
ChiralAM2AMx	IA +			
ChiralAM4AM4	IF + IE			
ChiralAM4AM5	IE + IG			
ChiralCE2CE4	IB + IC			
ChiralCE2CE5	IB + IK			
ChiralCE3CEx				
...				
ChiralAMCD-2	IA +	+ Astec Cyclobond I 2000 DMP		
ChiralAMCD-3	IF +			+ ChiralDex series LC columns
ChiralAMCD-5	IG +			+ ChiralDex series LC columns
CECD2	IB +	+ Astec Cyclobond I 2000 DMP		
CECD5	IK+			+ ChiralDex series LC columns
...				

Remark : PS: x =1-13, This table only reflects the compatibility between ChiralTek and other brands. The compatibility between other brands may not be correct in this table. Currently, ChiralTek has the most types of chiral products. The blank represents that other brands do not have this type of chiral column. For subsequent updates, please consult our local distributor.



www.chiralktek-column.com